

Recombinant SARS-CoV-2 NP NTD

domain V2

Catalog # EPT233

Expression Host E.coli

DESCRIPTION Recombinant SARS-CoV-2 NP NTD domain is

produced by our E.coli expression system and the

target gene encoding Gly44-Ser180 is expressed with

a 6His tag at the N-terminus.

Accession QHD43423.2

Synonyms 2019-nCoV coronavirus NP Protein; 2019-nCoV np

Protein; 2019-nCoV novel coronavirus Nucleoprotein

Protein

Mol Mass 18.7kDa

AP Mol Mass 18kDa, reducing conditions

Purity Greater than 95% as determined by reducing

SDS-PAGE.

Endotoxin

FORMULATION Supplied as a 0.2 µm filtered solution of PBS, pH 7.4

RECONSTITUTION



www.elkbiotech.com



SHIPPING

STORAGE

BACKGROUND

The product is shipped on dry ice pack.Upon receipt, store it immediately at the temperature listed below.

Reconstituted protein solution should be stored at ≤ -20°C.

Coronavirus N protein is required for coronavirus RNA synthesis, and has RNA chaperone activity that may be involved in template switch. Nucleocapsid protein is a most abundant protein of coronavirus. N protein packages the positive strand viral genome RNA into a helical ribonucleocapsid (RNP) and plays fundamental role during virion assembly through its interactions with the viral genome and membrane protein M. Plays an important role in enhancing the efficiency of subgenomic viral RNA transcription as well as viral replication. Because of the conservation of N protein sequence and its strong immunogenicity, the N protein of coronavirus is chosen as a diagnostic tool.





