



ELK Biotechnology

FITC-conjugated Goat Anti-Mouse IgG

Catalog NO. : SA008

For research use only.

Overview

Product name	FITC-conjugated Goat Anti-Mouse IgG
Recommended dilutions	1:20-1:100 for most applications NOTE: Optimal dilutions should be determined by the end user.
Size	100µL
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Note	This product is for research use only, not for diagnostic or therapeutic use.

Fluorophore FITC

Absorption Peak 492nm

Emission Peak 520nm